

E-tongue enables similarity / dissimilarity identification of wider range of products


Sensory panel limitations

- Training
- Availability
- Fatigue
- Food safe product



E-tongue

- ✓ More samples in a session
- ✓ No need for food grade product
- ✓ Close to continuous operation






Electrode array

BUT

Human sensory system is complex




16 electrodes coated with selective bio-polymers as proxy to physiological sensation \longrightarrow 15 voltage signals



Sensor

Reference solution

Test solutions



V7-V9


V8-V10 V9-V11


V10-V12 V11-V13 V12-V14 V13-V15

V14-V16 V15-V17 V16-V18

V17-V18

15 time series each transformed into 6 defined features






90 variables with many correlations



Random Forests variable importance shows number of features can be reduced


- Gini index distribution on 100 random forests
- 4 batches recorded: model trained on 3 (leave one out)
- Arbitrary choice of a cut-off limit at 1 in at least one batch
- F1, F2 and F5 for each of 15 series sufficient (blue vs black box plots)




Features

F11F21F5 only

Angle Classifier can be trained on a single category Is product under test IN or OUT of the category?






Positive: Labelled as OUT Sensitivity > Specificity



Dimension reduction and angle classifier enables good prediction accuracy with only data from a reference population



100% sensitivity with > 90% specificity

