

Adding the missing audit trail to R

Magnus Mengelbier
Managing Director

Topics

- Audit Trails
- Alternative approaches
- Tracking Changes
- Generating the Audit Trail
- Conclusion

Audit Trail

- Broader scope for and use of R with Life Sciences
 - Scope and use brings expectations on compliance and regulatory
 - Concept of Audit Trail established since the 1990's
 - Regulations
 - Directives
 - Guidance & Guidelines
 - R as a file-based application
- Independent record of who did what when
- ... Audit Trail tracks files that are created, updated and deleted

Alternatives

- Redefine R functions
 - Add audit trail methods to functions that write and delete files
 - Excludes any writes and deletions performed using external routines, e.g. c/c++, Python, Java, etc
- File system auditing and monitoring utilities
 - Records actions on the files system
 - Conflicts with IT security processes
 - Difficult to isolate R-related processes
- Repositories
 - Solution external to the R environment
 - Introduces new and complex compliance questions, challenges and issues
- Super log
 - Extend with logging functions to capture additional audit details

Tracking Changes

```
before <- list.files( ".", recursive = TRUE, include.dirs = FALSE)
```

```
#### <---- execute R program here
```

```
after <- list.files( ".", recursive = TRUE, include.dirs = FALSE)
```

```
created <- after[ ! after %in% before ]
```

```
deleted <- before[ ! before %in% after ]
```

	Prior	After
Created		x
Updated	x	x
Deleted	x	


Tracking Changes: Adding Updates

```
sapply( list.files( ".", recursive = TRUE, include.dirs = FALSE),  
        function(x) {  
            digest::digest( x, file = TRUE, algo = "sha256")  
        }, USE.NAMES = TRUE )
```

	Prior	After
Created		x
Updated	{digest}	{digest}
Deleted	x	

- Message Digest / Hash / Checksum is a derived value on the file content
- Different values indicates content is different
- Updated file = Different digest value before and after program execution

Audit Method



Executing R Programs: A Simple Approach Using Annotations

Annotate the program

```
# --- rexec controls
# rexec-clone /path/to/project/analysis
# rexec-inputs /path/to/project1/analysis/data adsl adae
# rexec-inputs /path/to/project2/analysis/data adsl adae
# rexec-outputs /path/to/project/analysis/data
```

Run the program

```
cxlib_rexec( "program.R" )
```


In Summary

- Audit Trails defined in industry regulations, directives and guidelines
- Different approaches inherits different complexities and challenges
- Audit Method
 - Plain R
 - Minimal set of additional packages
 - Containers and traditional server-based environment architectures
 - Any Operating System
 - Compatible with interactive R environments
- Extendable to other languages

Magnus Mengelbier

magnus.mengelbier@limellogic.com

Limellogic AB

Jungmansgatan 12
SE-211 11 Malmö
Sweden

www.limellogic.com
github.com/limellogic